

YEP Medium Version 2

Daniel Vaultot

Abstract

Use for *Schizochytrium*

Citation: Daniel Vaultot YEP Medium. [protocols.io](https://www.protocols.io)

<https://www.protocols.io/view/yep-medium-s2zegf6>

Guidelines



Before start

Necessary equipment

- Autoclave
- Laminar flow cabinet
- Stainless Steel Filter Holder
- Peristaltic pump


Solutions

- Seawater
- Nutriments (see protocol)

Plasticware and filters

- Polycarbonate bottle (Nalgene) : 1L
- Pipette
- Glass fibre prefilters (Millipore, AP1507500)
- Filters 0,22µm GSWP (Millipore, GSWP09000)
- Stericup® Filter Unit (Millipore, SCGVU10RE)

Materials

 Difco Bacto peptone [view by Fisher Scientific](#)

✓ D-glucose anhydrous by [Contributed by users](#)



Bacto Yeast Extract view by Becton-Dickinson

Protocol

Step 1.

Filter 1L of aged seawater (at least two months) on prefilter and 0.2 μm filter

Step 2.

Heat seawater during 20min at 100°C

 DURATION

00:20:00 Additional info:

Step 3.

To 300 mL of seawater, add :

- 1 g Difco bacto-peptone
 - 1 g Bacto yeast extract
 - 10 g D-glucose anhydrous
-

Step 4.

Complete final volume to 500mL of seawater

Step 5.

Autoclave the medium

Step 6.

Under laminar flow hood, filter the medium on 0.2 microns Stericup
