

Preparation Rice medium

1. Filter 1L of old seawater of at least two months on prefilter and filter 0,2 microns
2. Heat seawater during 20min at 100°C
3. To seawater, add 40 rice grains
4. Autoclave 3 times the medium (let temperature decrease between each autoclave)

The RCC team



Station Biologique de Roscoff,
CNRS/Université P. & M. Curie,
Place G. Teissier,
29680 Roscoff, France

*Apparatus : Autoclave, Laminar
flow cabinet, Stainless Steel
Filter Holders*

Solutions : Seawater, Rice grains

*Plasticware and filters : Bottle in
polycarbonate Nalgene, Glass
fibre prefilters (Millipore,
AP1507500), Filters 0,22µm
GSWP (Millipore, GSWP09000)*

ASSEMBLE

ASSOCIATION OF EUROPEAN MARINE BIOLOGICAL LABORATORIES