

Bold SW 1 NV modified

Adapted from [CCAP](#) recipe for use in seawater

Volume for 500 ml	Component	Concentration Stock Solutions		
5 ml	CaCl ₂	CaCl ₂ .2H ₂ O	0.5 gr	200 ml
5 ml	NaNO ₃	NaNO ₃	16 gr	200 ml
5 ml	MgSO ₄	MgSO ₄ .7H ₂ O	10 gr	200 ml
10 ml	PO ₄	From PCRS11 media		
5 ml	NH ₄ Cl (K)	From K media		
3 ml	K trace metal	From K media		
50 µl	Vitamines			
467 ml	SW sterile			